









Swab Guidelines

	<p>Amies Swab</p> <ul style="list-style-type: none">• Wound/Bacterial Culture: 10015• Throat Culture: 11234• Group B Strep: 43020
	<p>VCM Swab</p> <ul style="list-style-type: none">• Flu by PCR: 11404• RSV by PCR: 114042• Respiratory Viral Panel by PCR: 11204• Respiratory Bacterial Panel by PCR: 60010• Herpes simplex virus culture: 2649• Covid-19 by PCR: 74841 (Saline is an acceptable substitute)
	<p>BBL Culture swab (no gel or liquid media)</p> <ul style="list-style-type: none">• MRSA by PCR: 29212• Throat Culture: 11234
	<p>BD ESwab or Copan eSwab</p> <ul style="list-style-type: none">• Strep by PCR: 1010
	<p>BD Max UVE Specimen Collection Kit</p> <ul style="list-style-type: none">• CT/NG by PCR: 41378• Vaginal panel by PCR: 41390• MRSA by PCR: 29212 (Dry swab only, don't place the swab in UVE tube)
	<p>Aptima swab</p> <ul style="list-style-type: none">• HSV mRNA by PCR: 90570